

MOUTH PURE Oral Care Sponges

Features

- * Shape : Easy to take the dirt (Cherry blossom shape sponge)
- * Color : Easy to understand bleeding
- * Stick : Enables clean up every corner of the mouth by its long stick
- * Package : Individually packed, hygienically clean

How to Use

- * Remove from the bag, and moistened with a small amount of water in the sponge part, firmly squeezed.
- * Wipe the dirt of the oral cavity (surface of the teeth, gums, palate part, tongue etc.)
- * Enables to remove dirt effectively by rubbing the front from the back, turning round and round.

Size

- * Small (15mm dia.) Paper stem / Plastic stem
- * Medium (17mm dia.) Paper stem / Plastic stem
- * Large (20mm dia.) Paper stem / Plastic stem

